

STEAMJET®

**Class B**  
AUTOCLAVE

15  
min



Total process duration for 134°C process setting is just 15 minutes. Ultra short sterilization time allows preparing tools when they are needed.

**2 YEARS WARRANTY** **24<sup>7</sup> SERVICE**



The large touchscreen facilitates intuitive use. In addition to basic functions, it is possible to set up your own sterilization cycle.

The process is presented in a graphic form giving information on remaining time, temperature and pressure.



STEAMJET®



**2 YEARS WARRANTY** **24<sup>7</sup> SERVICE**


All STEAMJET autoclaves are covered by a full three year warranty available on-site. ENBIO Service performs warranty claims 24 hours a day, 7 days a week.

Program	121°C	134°C
Tools	wrapped and unwrapped	
Sterilization time	15 min	4 min
Drying time	5 min	4 min
<b>Total process duration</b>	<b>25 min</b>	<b>15 min</b>

**Class B**  
AUTOCLAVE

TECHNICAL DATA

Power supply	230V/50 Hz
Power output	3500 W
Maximum current consumption	15 A
Operating pressure	2,1 bar
Maximum pressure	2,3 bar
Maximum process temperature	137°C
Volume of the process chamber	2,7 dm <sup>3</sup>
Process chamber dimensions (LxWxH)	292x192x45 mm
External dimensions (LxWxH)	561x252x162 mm
Weight	15 kg
IP Code	IP20
Process data logging	USB Drive



Enbio Technology Sp. z o.o.  
T. +48 58 731 96 31  
info@enbiogroup.eu  
www.enbiogroup.eu

enbio®

The partnership is registered in the register of the entrepreneurs, kept by the District Court for Gdańsk-Północ, Commercial Division 8 of the National Court Register, company registration number 214555, National Official Register of Business Entities Number 191317763, Tax Identification Number 958-09-82-035. The capital paid in full is 16,763,087.29 PLN (which is the equivalent of EUR 4,039,298.14, with the exchange rate 4.15 EUR/PLN), including share capital of 887,500 PLN.

STEAMJET®

**Class B**  
AUTOCLAVE

Autoclave for instruments



www.enbiogroup.eu

enbio®

STEAMJET®

Compact

Ultraquiet

Rapid

Intuitive

## Sterilization at hand

Compact design, small footprint and extremely quiet operation allow you to use autoclave directly in the office. Instant sterilization process makes each instrument available just in time.

